





8	It is possible to mount only the <b>FP-05</b> chest belt equipped with buckles, fastened circumferentially. The chest belt equipped with the mounting triangles cannot be mounted.
9	In case of the <b>JORDI™</b> chair, it is possible to mount the <b>FP-21</b> belts in the narrow ( <b>FP-21/W</b> ) as well as the wide ( <b>FP-21/S</b> ) version, provided that the seat is not equipped with a footrest. In case of the <b>JORDI™</b> chair equipped with a footrest or in case of the <b>JORDI HOME™</b> stroller (which is equipped with a footrest as standard), it is suggested to buy a narrow version of the belt.
10	To mount the vest, you must buy the "female" buckle member (mounting element <b>FP-EM-8</b> ) and mount it on the vest as instructed. <b>ATTENTION!</b> It is possible to mount only the vest with the rear adjustment <b>T</b> – caregiver side. After mounting the vest, the adjustment is no longer available. It is impossible to mount the vest with the front adjustment <b>P</b> – user side.
11	Mount the belts with the use of the wood screw ( <b>FP-EM-6</b> ) in the designated place. You need to buy <b>two pieces</b> of the element <b>FP-EM-6</b> .
12	Mount the belts in the designated place with the use of the mounting element: - For <b>KIDOO™</b> size 1, 2, 3 – <b>FP-EM-4-M6x16</b> - For <b>KIDOO™</b> size 4, 5 – <b>FP-EM-4-M6x20</b> You need to buy <b>two pieces</b> of the mounting element.
13	Mount the belts in the designated place with the use of the mounting element: - For <b>KIDOO™</b> size 1, 2, 3 – <b>FP-EM-4-M6x16</b> - For <b>KIDOO™</b> size 4, 5 – <b>FP-EM-4-M6x20</b> You need to buy <b>four pieces</b> of the mounting element.
14	Ankle stabilizers are available with the <b>KIDOO HOME™</b> and <b>JORDI HOME™</b> strollers, and also with the <b>KIDOO™</b> and <b>JORDI™</b> chairs equipped with footrests (for <b>KIDOO™</b> it is <b>KDO_616</b> or <b>KDO_618</b> ; for <b>JORDI™</b> it is <b>JRI_619</b> ). The size of the ankle stabilizers ( <b>FP-22</b> ) should be selected based on the measurement of the leg above the ankle, according to the size table on the website <a href="http://flex-point.pl/en/">flex-point.pl/en/</a> .
15	It is recommended to use the extender for universal stabilizing belt ( <b>FP-A-1</b> ).
16	The size of the ankle stabilizers ( <b>FP-22</b> ) should be selected based on the measurement of the leg above the ankle, according to the size table on the website <a href="http://flex-point.pl/en/">flex-point.pl/en/</a> .

**ATTENTION!** The sizes indicated in the table are only a suggestion.

Before buying the **FlexPoint®** belts, check the size table on the website [flex-point.pl/en/](http://flex-point.pl/en/) and select the appropriate size of the **FlexPoint®** product for the user.

**ATTENTION!** The detailed instructions for mounting the **FlexPoint®** products can be found in the user's manual of a given stroller or chair.

The user manuals are available on the website [en.akces-med.com](http://en.akces-med.com) in the given product tab.